

ADHESIVE MORTAR

**PRODUCT NAME**

CLINK' BASE

STANDARD

UNE - EN 12004

COMPOSITION

Grey/white cement, aggregates and additives.

DESCRIPTION

CLINK' BASE is thin-coat adhesive mortar for use with medium/high porosity tiles on both floors and interior walls. Available in grey or white





ADHESIVE MORTAR

CLINK' BASE

METHOD OF APPLICATION

PREPARATION:

- Ensure the surfaces to be tiled are in good condition and clean, that there are no traces of grease or dust, and that they are mechanically resistant.
- The plane should differ by no more than 5mm over a distance of 2m.

PRODUCT PREPARATION:

- Sprinkle the required amount of CLINK' BASE into clean water.
- Mix by hand or with an electric mixer (on a low setting), until you have a smooth, even mixture with no lumps.
- Leave to stand for around 2 minutes to let it hydrate correctly

APPLICATION:

- Use a trowel to apply a thick layer of CLINK' BASE to the prepared surface and then spread it out evenly with a toothed trowel.
- Place the tiles on top and press until the grooves in the mortar have been squashed.



ADVANTAGES

High yield

Suitable for all types of traditional tiling

High yield

Excellent adhesion

TECHNICAL DATA

WATER / PRODUCT RATIO:	6 - 8 l/sack
MIXTURE STANDING TIME:	2 min
CORRECTION TIME:	15 min
USE WITHIN:	3H
WALK ON AFTER:	24H
YIELD:	4 - 5 kg/m ²
MIXING SPEED:	500 rpm
APPLICATION TEMPERATURE:	5 °C - 40 °C

PACKAGING



20 KG/SACK



64 SACKS/PALLET



25 TN



STORAGE AND PRESERVATION

- Store the product on a flat surface and never outdoors. Keep dry and out of direct sunlight.
- Use within 12 months of the date of manufacture.

LEGAL NOTICE: This information, in particular the recommendations relating to product application and the technical data, is based on our experience and current knowledge as well as the suggested uses and common applications of our products, and is subject to the final conditions of the site or the application of the products. GYPROD is not liable for the use of the products by third parties, or for the results obtained.
The GYPROD quality guarantee is limited exclusively to the product supplied.