

FLEXIBLE JOINT MORTAR



PRODUCT NAME	CLINK' JOINT FLEX
STANDARD	UNE – EN 13888
COMPOSITION	Grey/white cement, sand, resin and various additives.
DESCRIPTION	CLINK' JOINT FLEX is a flexible mortar for tiling joints on indoor and outdoor walls and floors. It is particularly suitable for grouting tiles on deformable surfaces (wooden floors, heated floors, etc.) Available in grey or white.



FLEXIBLE JOINT MORTAR

CLINK' JOINT FLEX

METHOD OF APPLICATION

PREPARATION:

- The joints must be clean and dust-free.
- The tiles must be laid at least 24 hours before jointing is carried out.

PRODUCT PREPARATION:

- Sprinkle CLINK' JOINT FLEX into clean water. Mix by hand or with an electric mixer (on a low setting), until you have a smooth, even mixture with no lumps.
- Leave the mixture to stand for 10 minutes.

APPLICATION:

- Use a rubber trowel to fill the joints by spreading the product diagonally and pressing it into them.
- Once the mortar has begun to harden, it can be cleaned off with a damp sponge.
- The final cleaning should be done when the product has hardened completely by rinsing the surface with clean water or wiping it with a damp cloth.

ADVANTAGES



High yield

High deformability

Excellent adhesion

Maximum safety

TECHNICAL DATA

WATER / PRODUCT RATIO:	5 l/sack
MIXTURE STANDING TIME:	10 min
MIXING SPEED:	500 rpm
DRYING TIME:	24H
APPLICATION TEMPERATURE:	5 °C - 40 °C

PACKAGING



10 KG/SACK



120 SACKS/PALLET



25 TN



STORAGE AND PRESERVATION

- Store the product on a flat surface and never outdoors. Keep dry and out of direct sunlight.
- Use within 12 months of the date of manufacture.

LEGAL NOTICE: This information, in particular the recommendations relating to product application and the technical data, is based on our experience and current knowledge as well as the suggested uses and common applications of our products, and is subject to the final conditions of the site or the application of the products. GYPROD is not liable for the use of the products by third parties, or for the results obtained. The GYPROD quality guarantee is limited exclusively to the product supplied.